



Disinfectants for TUFFease Products

Listed below are common medical facility cleaning/disinfectant solutions. Active ingredients and acceptability are noted. The spray-wipe-rinse procedures are assumed with a thorough clear water rinse as a final step. Wash temperatures up to 103F/40C are acceptable. Do not soak. Gas sterilization, autoclaving, machine washing and drying aren't recommended.

Brand Name	Active Ingredient	Acceptable
2030 Bio-Cidol [®] Neutral	Dymethyl Ammonium Chloride	Yes
Accel	Accelerated Hydrogen Peroxide	Yes
Antifect AF	Dymethyl Ammonium Chloride	Yes
Betadine	Povidone-Iodine	Yes
Bruphene 256	Amylphenol	Yes
Buraton 10F	Dialdehyde	Yes
Clorox Bleach	Sodium Hypochlorite, 6%	Yes
Defend PS-10	Dymethyl Ammonium Chloride	Yes
Dispatch	Sodium Hypochlorite, 1%	Yes
Eco Tru	Chloroxylenol	Yes
Endbac	Dymethyl Ammonium Chloride	Yes
Hebiclens Liquid	Chlorhexidine Gluconate	Yes
PDI Sani-Cloth [®] Plus	Dymethyl Ammonium Chloride	Yes
Phenex N. S.	Amylphenol	Yes
Quat-64 Broad	Dymethyl Ammonium Chloride	Yes
Spartan TB-Cide Quat	Quaternary Ammonium	Yes
Terralin	Dymethyl Ammonium Chloride	Yes
Vesphene LpH se	Phenylphenol	Yes
Virex TB	Dymethyl Ammonium Chloride	Yes
CaviCide	Isopropyl Alcohol	No
Endozime AW	Isopropyl Alcohol	No
Mikrozyd Liquid	Isopropyl Alcohol	No

Never use any type of alcohol solutions to clean Tuffease[®] pads

This chart is for general reference and a guideline only. Check with the solution manufacturer for specific application regarding viruses targeted and effective exposure times required.