

# *imp* De Mayo Knee Positioner® Set Up Instructions

1. Replace OR Table Pads with green *imp* pads




If using X-Ray cassette, it is not necessary to remove the X-ray cassette plates.



2. Position patient with gluteal fold at separation of pads



3. During prep of surgical leg, remove single OR Table Pad



4. Create a 3" - 5" fold in final drape under the buttocks



5. Place De Mayo Knee Positioner® in the well to achieve maximum flexion





6. Check flexion for final approval of freedom of the drape to allow full range of motion when the Positioner is locked on drapes



7. Position clamp directly under gliding pins and lock to the OR Table Side Rail with the Single Lever Clamp

*Single Lever Clamp Set Up*




1. Position the clamp holding with left hand and with the lever in "unlock" position



2. Rotate lever to "lock" position



3. Squeeze jaws tight on side rail, and tighten knurled knob all the way





### Single Lever Clamp Set Up Continued



4. While holding clamp, rotate lever to "unlock"




5. Turn knurled knob 1/4 to 1/2 rotation



6. Rotate lever to "lock" position requiring only 2 finger pressure to achieve locking. Works like a vice grip and will "snap" into place



7. Slide pins into clamp holes





After seating pins into holes of the clamp,



tighten knurled knobs



Single Lever Clamp Set Up is complete

8. Protect your patient from potential Pressure Ulcers. Insert the *imp* Patient Protective Pad into the sterile boot



**WARNING**  
Liability Voided if *imp* Patient Protective Pads are not used.

9. Wrap cohesive bandage around the foot starting with a minimum of six (6) foot wraps, tear, and finish the wrap above the bracket around the calf.

*Avoid wrapping cohesive bandage over the bracket on the back of the boot.*

