

*imp* **UNIVERSAL LATERAL  
POSITIONER™ SYSTEM**

*We've Made Adjustments  
So You Won't Have To*

*PROVIDING MORE POSITIONING  
OPTIONS TO ACCOMMODATE A  
RANGE OF PATIENT ANATOMIES*

US Patent No. 6,003,176

**IMP...THE OPERATIVE WORD  
IN PATIENT POSITIONING**



# UNIVERSAL LATERAL POSITIONER™ SYS

## MORE POSITIONS FOR BETTER POSITIONING

The Universal Lateral Positioner™ is ideal for any surgery requiring the patient to be fixed in a stable lateral decubitus position. It provides a reproducible point of reference for primary and revision arthroplasty. The positioner is available in two sizes, a 22" base plate model, and a 25" base plate model. Because obesity poses its own set of problems when performing surgery, the Large Universal Lateral Positioner (25"), Cat. #1092L, has an additional 3" added to the base plate and 1-1/2" to the anterior arms to accommodate fatty tissue.



Redesign of the anterior support arm provides easy access to the locking/adjusting knobs.



Removable and independently adjustable anterior arms allow more accurate contact with the iliac crests.



Inferior arm is removable to accommodate all patient anatomies, while still providing effective stabilization in the lateral decubitus position.



Superior arm on the 25" model adjusts horizontally allowing for patient girth or positioning in a plane other than vertical.



Patented Hyperflexion Plate, Cat. #298, allows up to 120° of flexion during test for range of motion.

# SYSTEM



Positioner without  
Hyperflexion Plate

Low profile side rail clamp  
included in System

# imp UNIVERSAL LATERAL POSITIONER™ SYSTEM



Contour Lumbar Pad, Cat. #1002-P, conforms to patient anatomy improving alignment and pressure distribution.



Latex Free Foam Pad Set



Dual Compression Foam Pads for anterior support arms are unsurpassed in preventing skin irritation.



Storage Case keeps components organized and accessible for OR staff while protecting capital equipment investment. Sturdy casters provide easy handling.

## Ordering Information:

### Description

Universal Lateral Positioner System

Includes 1 ea. Universal Lateral Positioner – 22” (1092), Hyperflexion Plate (298), OR Table Rail Clamp (405-RC), and Storage Case (1092-Case)

Universal Lateral Positioner System – Large

Includes 1 ea. Universal Lateral Positioner – 25” (1092-L), Hyperflexion Plate (298), OR Table Rail Clamp (405-RC), And Storage Case (1092-Case)

Universal Lateral Positioner – 22”

Universal Lateral Positioner – 25”

Hyperflexion Plate

OR Table Rail Clamp

Storage Case w/Casters

Contoured Foam pads (6 sets per/case)

Disposable Foam Pads (6 sets per/case)

Latex Free Gel Pads

### Cat. No.

1092-SYS

1092L-SYS

1092

1092-L

298

405-RC

1092-CASE


1002-P

1092-P

1092-GP



Distributed By:



**imp** Innovative  
Medical Products, Inc.

*A GLOBAL LEADER IN PATIENT POSITIONING*

87 Spring Lane, Plainville, CT 06062

PH: 860-793-0391 FAX: 860-793-8975

Toll Free: PH: 800-467-4944 FAX: 888-229-1452

[www.impmedical.com](http://www.impmedical.com)

US Patent No. 6,003,176

Universal Lateral Positioner is a trademark of Innovative Medical Products, Inc.

5290 10M ©2005 Innovative Medical Products, Inc. All rights reserved. 12/05