

imp® Humbles LapWrap® Positioning Pad


PATIENT SAFETY DEVICE IDEAL FOR ROBOTIC LAPAROSCOPY



"For me, securing a patient's arms by their sides has always been a real pain. Since using the LapWrap®, it has become a piece of cake. I like the ease of use, the protection it affords and the fact that I can access the patient's arms during the case without interrupting surgery. When I have the need to tuck the patient's arm, one of the first things I do is make sure the LapWrap® is on the OR table. My only complaint is that it was not available when I first went into practice!"

Thom L. Holt, CRNA
Conway, South Carolina

- Allows easy access for IV tubes and Leads for anesthesia
- Secures patients' arms by their sides; helps prevent shoulder injury by *NOT* allowing a patient's arm to fall from the side of the OR table
- For laparoscopic, head/neck and neurosurgery
- Protects neurological structures by preventing hyperextension of the elbow
- Universal Design Fits All Size Patients
- Helps retain the patient's body heat



For unusual table positions, place center tongue between OR table and side rail.



Optional Technique: Place TUFFease® pressure redistribution pad on OR table for extra patient protection.



Cat.#: 904-G
10/Case



U.S. Patent No. 8,001,635

Always use appropriate positioning safety precautions

Latex Free



LapWrap® is a registered trademark of Innovative Medical Products, Inc.
LapWrap® - U.S. PATENT NO. 8,001,635
TUFFease® is a registered trademark of Innovative Medical Products, Inc.
© 2012 INNOVATIVE MEDICAL PRODUCTS, INC. ALL RIGHTS RESERVED

imp® Innovative
Medical Products, Inc.
A GLOBAL LEADER IN PATIENT POSITIONING

87 Spring Lane, PO Box 8028
Plainville, CT 06062
PH: (860) 793-0391 FAX: (866) 459-1805
Toll Free PH: 1-800-467-4944 FAX: 1-888-229-1452
www.impmedical.com