

imp De Mayo Knee Positioner®



" Like having another pair of hands in the O.R."

***Single Locking Lever Allows for Precise Control
of Flexion, Extension, Tilt and Rotation***

impTHE OPERATIVE WORD
IN PATIENT POSITIONING

De Mayo Knee Positioner®

The De Mayo Knee Positioner® reduces surgeon fatigue. The unique **single locking lever** allows for precise release of the boot for testing of ligament in the sterile field on top of the sterile drapes.




and assists in the stabilization of a patient's leg during surgery. control and its ease of use saves surgical time. The same lever stability and ROM testing. The positioner is easily set up




The **single locking lever** controls flexion, extension and boot release; as well as finite internal/external rotation allowing superior access to the lateral and medial compartments.

Single Lever Clamp securely attaches positioner to O.R. table siderail over the sterile drape without perforation. (Cat. #903)



External Rotation



Internal Rotation



Flexion



Extension

imp® De Mayo Knee Positioner® Accessories



5 Section OR Table Pads
(Cat. No. 398-LG Kit)



(Cat. No. 294-P-10)

Single use, latex-free sterile foam pad and cohesive wrap protect the patient from pressure ulcers, abrasions and possible neurological impairment while securing the foot into the boot.

De Mayo Universal Distractor™



(Cat. No. 907)




(Cat. No. 907-P)

Sterilization Cases Kit



Accessory Case



Base Case

(Cat. No. 706 Kit)

Ordering Information:

Description

De Mayo Knee Positioner® 30" Base Plate, Aluminum Boot & Single Lever Clamp

De Mayo Knee Positioner® 25" Base Plate, Aluminum Boot & Single Lever Clamp

Accessories:

De Mayo Knee Positioner® Aluminum Boot with Distractor Block

De Mayo Knee Positioner® Composite Boot with Distractor Block

De Mayo Universal Distractor™

Center Clamp Threaded

Single Lever Clamp

5 Piece Set O.R. Table Pads for use with De Mayo Knee Positioner®

(includes 3 piece O.R. Table Pads & 2 piece Armboard Pads)

Sterilization Cases Kit

(includes Accessory Case & Base Case)

Patient Protective Pads

De Mayo Knee Positioner® Sterile Foam Pads & Cohesive Wrap 10/Case

De Mayo Universal Distractor™ Sterile Pressure Protective Pad Set 10/Case

(includes Distractor Pad, Foam Pad & Cohesive Wrap)

Cat. No.

803 Kit

803-S Kit

803-ABD

803-CBD

907

302

903

398-LG Kit

706 Kit

294-P-10

907-P



De Mayo Knee Positioner® is a registered trademark of Innovative Medical Products, Inc.
De Mayo Knee Positioner® U.S. Patent # 7,380,299
De Mayo Universal Distractor™ is a trademark of Innovative Medical Products, Inc.
De Mayo Universal Distractor™ U.S. Patent # 8,048,082
© April 2012 Innovative Medical Products, Inc. ALL RIGHTS RESERVED

imp® Innovative
Medical Products, Inc.

87 Spring Lane, PO Box 8028
Plainville, CT 06062

PH: (860) 793-0391 FAX: (866) 459-1805

Toll Free PH: 1-800-467-4944

www.impmedical.com