



imp

PORTA-VIEW[®]

Portable X-Ray View Box



Innovative Medical Product's Porta-View™ conveniently displays images from X-ray or MRI scans at the operative site. The maneuverability of the view box allows the surgeon to view films without leaving the sterile field. The Porta-View™ eliminates the re-adjustment of magnifying lenses or loupes.

The tilt and height are adjustable for easy viewing while sitting or standing and viewing angle adjusts from 0° to 90°. An optional sterile transparent cover can be utilized to trace the x-ray, design an osteotomy, or decide on hardware placement.

PORTA-VIEW™

- 28" x 17" viewing area
- Sterile transparent cover (optional)
- Adjusts viewing angle 0° to 90°
- Viewing height from 36" to 60"
- Sturdy, maneuverable stand



Surgeons in the fields of Maxiofacial Oral, E.N.T., Neurology, Orthopedic, (i.e. foot and ankle, hand and wrist), Pediatrics, Radiology and other disciplines can benefit from its use.

Gravity roller grip retainers and 25 ft. cord with a hospital grade plug allows versatility and efficiency.

© Copyright 1998 Innovative Medical Products, Inc. All rights reserved.

Innovative design of the Porta-View™ by G. David Rutland, MD. based on years of surgical experience.

imp[®] **innovative**
MEDICAL PRODUCTS, INC.

INNOVATIVE MEDICAL PRODUCTS, Inc.
87 Spring Lane, Plainville, CT 06062
Tel: 800.467.4944 860.793.0391
Fax: 800.793.8975
Toll Free Fax: 888.229.1452
Email: info@impmedical.com
www.impmedical.com

Distributed By:

